

Bath Transfer Systems Measurement Guide

Bath/ commode chairs and slider systems are known for their ability to fit in cramped spaces others simply can't. While we have yet to find a bathroom we can't accommodate, a few simple measurements will insure the best fit of your chair. Follow the descriptions below, and enter the measurements. Round up or down to 1/2-inch increments, and use the decimal system (for example, 10.5 for 10 and-half inches).

A Overall length of bathtub:

B Length of bathtub floor:

C Width of bathtub floor (flat portion only - do not include any rounded area)

D Overall width of bathtub:

E Distance between edge of toilet bowl and bathtub:

F Widest portion of toilet bowl (including "wedge" shape up to face tank)

G Distance between edge of toilet bowl to closest objects/wall:


H Distance from front edge of toilet bowl to closest objects/wall:

K Distance from center of toilet bowl to center of bathtub drain/spout:

L Distance from top edge of bathtub to floor of bathtub:

M Distance from top edge of bathtub to bathroom floor:

N Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)



Note: Please verify the measurements of the product, to ensure proper fitting for the client and adequate space availability in the user environment. medicashop will not be responsible for incorrect measurements/fitting and inadequate space availability in the user environment.

Need help with measurements, Please call us...

HOME HEALTH PAVILION

5027 Route 9 W Newburgh,
NY 12550
Phone: (845) 569 1250
Fax: (845) 569 1291

87 Danbury Road, Unit I New Milford,
CT 06776
Phone: (860) 210 1313
Fax: (860) 354 1123

Email : sales@homehealthpavilion.com