


Seating Measurement Guide



A Seat Pan to Top of Shoulder: _____

B Chest Depth: _____

C Chest Width: _____

D Seat Depth: _____
(Measurement from behind the knee)

E Seat Pan to Top of Head: _____

F Seat Pan to Elbow: _____

G Hip/Seat Width: _____

H Behind Knee to Heel: _____
(Lower leg length, seat pan to footrest)

I Seat to Axilla/Armpit: _____

Age: _____

User Weight (kg/pounds) _____

User Height _____

Note: Please verify the measurements of the product, to ensure proper fitting for the client and adequate space availability in the user environment. medicaleshop will not be responsible for incorrect measurements/fitting and inadequate space availability in the user environment.

Need help with measurements, Please call us...

HOME HEALTH PAVILION